

β -Gal Staining of *Xenopus* Embryos

- 1) Fix embryos for 30 minutes to 1 hour (no more!) in MEMFA (while nutating).
- 2) Rinse 2 times with PBS + 2 mM MgCl₂.
- 3) Transfer to staining solution (max. volume).
- 4) Transfer to X-gal mix: add 100 λ of the 50 mg/ml stuff (it is a small brown bottle in the -20 °C freezer) per 5 ml of embryo sample in staining solution. Use a minimal volume (5 ml is more than enough for 1 scint. vial) – this stuff is expensive. To develop more rapidly and with a deeper purple color, add 5 λ of NBT (NBT stock as for *in situs*: 75 mg/ml in 70% dimethyl formamide / 30% water) to 10 ml of final developing mix.
- 5) Incubate at **37 °C** in the dark (wrapped in tin foil) until desirable stain develops (a few hours to 24 hours).
- 6) Rinse in PBS, and twice in Methanol. Samples can be stored here.
- 7) Clear in Benzyl Benzoate : Benzyl Alcohol (2:1) if desired. BB:BA is very toxic; it is not water mixable (clean glass with ethanol) and melts plastic, so use it in glass only. It will also leech the stain out of embryos (which become brittle) in a while. You should photograph your samples in a watch-glass.

Staining solution (kept wrapped in foil in fridge at 4 °C):

	<u>10 ml =</u>	<u>1 L =</u>
20 mM K ₃ Fe(CN) ₆	0.06 g	6 g
20 mM K ₄ Fe(CN) ₆ ·3H ₂ O	0.08 g	8 g
0.01% deoxycholic acid	0.001 g	0.1 g
0.02% NP-40	2 λ	200 λ
in PBS + 2 mM MgCl ₂		